
A		B		C		D		E		F	
---	--	---	--	---	--	---	--	---	--	---	--

REV.	ZONE	DESCRIPTION	DCO NO.	DCO BY	UPDATED
001	N/A	INITIAL RELEASE	2017-MKT-DCO-042	S.A.	AUG 2017
002	N/A	ADDED DISPOSAL STATEMENT	2018-MKT-DCO-231	S.A.	OCT 2018
003	N/A	UPDATED PER INTERNATIONAL LABELING REQUIREMENTS	2018-MKT-DCO-068	S.A.	OCT 2018

1
2
3
4

- NOTES:
1. TO BE PRINTED ON WHITE 20 LB. BOND (OR EQUIVALENT).
 2. IFU TO BE A TWO SIDED BOOKLET.
 3. SIZE TO BE 8.5 +/- .25 HIGH X 11.0 +/- .25 WIDE.
 4. FOR ARTWORK SEE: 3000314-21 UNIVERSAL PLUS AND UNIVERSAL PLUS DUAL.indd
 5. IFU TO BE LEGIBLE AND FREE OF ANY SMEARS OR TEARS.
 6. LINES SHOWN FOR REFERENCE ONLY.
 7. SHEETS TO BE FASTENED WITH A SINGLE STAPLE IN THE UPPER LEFT HAND CORNER.

THIS DRAWING WAS RE-CREATED IN ADOBE INDESIGN

MEGADYNE [®]		MEGADYNE MEDICAL PRODUCTS INC. 11506 SOUTH STATE STREET DRAPER, UT 84020-9453 USA (801) 576-9669	
TOLERANCES UNLESS NOTED OTHERWISE .X = ± .1 .XX = ± .02 .XXX = ± .005 ANGLES = ± 1 FRACTIONS = ± 1/16	ENGINEER ME	DATE 29-AUG-2017	TITLE IFU, UNIVERSAL PLUS UNIVERSAL DUAL PLUS
	DESIGNER S.A.	DATE 29-AUG-2017	SIZE B
			PART NUMBER 3000314-01
			REV. 003
DIMENSIONS INCHES	ASME Y14.5 2009	DO NOT SCALE PRINT	SCALE N/A SHEET 1 OF 27

FRONT COVER

DIMENSIONS AND SPECIFICATIONS IDENTIFIED AS CONTROLLING CHARACTERISTICS ARE DENOTED BY "◆".
THIS DOCUMENT IS THE CONFIDENTIAL PROPERTY OF MEGADYNE MEDICAL PRODUCTS, INC. AND INCLUDES PROPRIETARY INFORMATION. IT IS TO BE USED FOR NO PURPOSE OTHER THAN THE LIMITED INTENT FOR WHICH IT IS DELIVERED AND NOT TO BE USED BY OR DISCLOSED TO ANY OTHER INDIVIDUAL OR LEGAL ENTITY WITHOUT THE EXPRESS WRITTEN CONSENT OF A MEGADYNE REPRESENTATIVE.

1	A	B	C	D	E	F
2		<p style="text-align: center;">ETHICON Megadyne™</p> <p style="text-align: center;">Megadyne Medical Products Inc. 11506 South State Street Draper, Utah 84020 USA Phone: +1 (801) 576-9669</p> <p style="text-align: center;">Universal Plus REF 0847 Universal Dual Plus REF 0848 Patient Return Electrode with Polymer "English"</p> <p style="text-align: center;">INDICATIONS FOR USE</p> <p>This device is designed to be used whenever monopolar electrosurgery is indicated. The intended use of this device is to conduct monopolar electrosurgical energy from target tissue of a patient back to one or two electrosurgical units (ESU), or generators.</p> <p>Electrosurgical use is restricted to use with isolated monopolar electrosurgical generators. The device is not intended for RF ablation.</p> <p style="text-align: center;">INSTRUCTIONS FOR USE</p> <p>The device is intended to be used as a 24 month reusable patient return electrode for all patients from .35 kg (0.8 lb) and up (no maximum restriction on weight and no restriction on age). The Universal Plus Patient Return Electrode has a single cord that connects to one generator. The Universal Dual Plus Patient Return Electrode has two cords that can connect to one or two generators.</p> <p>To Install the Universal Plus or Universal Dual Plus Patient Return Electrode:</p> <ol style="list-style-type: none"> 1. Inspect the Patient Return Electrode for damage to the outer skin, cable(s), or connector(s) (e.g. holes, tears, peeling, damaged insulation). Replace the Patient Return Electrode or cable(s) before using if electrical components or cable insulation are damaged or if the month of expiration is reached. Properly discard after use per institution's policy for disposal of biomedical waste. If the outer skin is damaged, repair before use. See repair instructions below. 2. Place the Patient Return Electrode on the operating surface, green side up. Do not place the Patient Return Electrode directly onto a metal surface. The Patient Return Electrode may be placed over positioning devices and frames. 3. Position the patient directly on the Patient Return Electrode to maximize patient contact area and for best skin shear relief. The Patient Return Electrode should not be used on top of the patient. Make sure the electrode remains smooth and prevent wrinkles or folding of electrode under patient. 4. If the Patient Return Electrode requires repositioning under the patient, please use care to not cause undue stress or stretching which may result in splitting or damage of the electrode. 5. Do not place excessive amounts of linens or other materials between the patient and the Patient Return Electrode. The use of excessive material between the patient and the electrode assembly may result in a diminished electrosurgical effect. A sheet and draw sheet may be placed over the Patient Return Electrode; however, this reduces the shear relieving characteristics. If a sheet is used, it should fit loosely over the electrode and be unstarched and free of wrinkles or folds. 6. To install the Patient Return Electrode, insert the correct cable connector (REF M2K-01, M2K01, M2K-02, M2K02, M2K-03, M2K03, M2K-04, M2K04, M2K-05, M2K05, M2K-06, M2K06, M2K-08, M2K08, M2K-09, or M2K09) for the generator being used. Ensure cable is not located under patient or where it could be damaged. 7. The Universal Dual Plus Patient Return Electrode may be used with two separate generators. When using only one generator with the Universal Dual Plus Patient Return Electrode, use the pre-attached cap(s) supplied with each cable to cover the exposed connector and secure the cord out of the way before starting electrosurgery. If the cap(s) is damaged or lost, call Megadyne customer service for a replacement. If the generator allows for a single or dual plate selection, select the single plate option. If the cable is installed correctly, the generator light will change from red to either no light or a green light depending on the model of generator. If the light does not change, DO NOT PROCEED. Contact your Megadyne customer service representative. 	<ol style="list-style-type: none"> 8. To install the Patient Return Electrode in a contact quality monitoring (CQM) dedicated generator, insert the blue cable connector (M2K-08, M2K08, M2K-09 or M2K09) end into the patient return electrode receptacle on the generator. The CQM light should go green, indicating that the electrode has been recognized by the generator and allowing for the use of electrosurgery. If the CQM light of the generator does not illuminate after connecting the electrode to the generator, DO NOT PROCEED. Contact your Megadyne customer service representative. Note: The Patient Return Electrode utilizes current limiting technology which has been determined to be safer than the contact quality monitoring system regarding pad site burns. <p>To disassemble the Patient Return Electrode:</p> <ol style="list-style-type: none"> 9. Turn the generator to the "OFF" position, or "STANDBY" position. 10. Disconnect the cable from the generator by unplugging it. Store the cables when not in use to prevent a tripping hazard or damage to the cables. 11. The electrode should be stored either on the operating room table or in an area away from heavy or sharp objects when not in use. <p>Cleaning and Care Instructions:</p> <ol style="list-style-type: none"> 12. Clean/disinfect the Patient Return Electrode, and cable with any of the following solutions: mild bleach solution (10:1); glutaraldehyde; o-phenylphenol; o-benzyl-p-chlorophenol; p-tertiary amylphenol. Thoroughly rinse the Patient Return Electrode with clear water to remove any residue from cleaning solutions. Failure to remove chemical build up may cause premature failure of the outer skin. <p>Alcohol and other strong disinfectants may cause staining or hardening of the electrode's outer skin. <u>Cleaning solutions with peroxide should not be used.</u></p> <ol style="list-style-type: none"> 13. Do not allow cleaning solutions or rinse water to come in contact with metal connectors. 14. Dry the return electrode prior to next use. <p>Transportation:</p> <ul style="list-style-type: none"> Do not pull, lift, or carry by the cord. This may cause electrical failure. Always handle the Patient Return Electrode with care. Carefully remove it from the table or other surgical equipment. The Patient Return Electrode should be loosely rolled or folded and cradled (like a baby) when carried. Carrying the electrode unsupported can stretch the outer skin and cause undue stress that may result in splitting. Lifting the product from the corners may stretch the outer skin, causing stretch marks or tears. <p>Repair:</p> <ul style="list-style-type: none"> A MEGADYNE Repair Kit and scissors can easily be used to repair damage to the electrode's outer film. Follow the instructions for use included with the repair kit for more information. If the Patient Return Electrode is not functioning as expected, a biomedical technician may measure the resistance of the cables and corner connection. If the resulting resistance is greater than 3 ohms end to end for all cables except the M2K-08, M2K08, M2K-09 or M2K09, please remove the electrode and cables from service and contact Megadyne customer service for Patient Return Electrode or cable replacement according to the terms of the warranty. For the M2K-08, M2K08, M2K-09 or M2K09 cables, a capacitive reading outside a range of 90 nF and 110 nF, end to end, remove the cable from service and contact your customer service representative for a replacement according to the terms of the warranty. 	<p>Heating and Cooling:</p> <ul style="list-style-type: none"> The patient may be placed on a pre-warmed Patient Return Electrode in order to minimize patient heat loss. The Patient Return Electrode may be pre-warmed by placing it under or over a patient warming device. The Patient Return Electrode is NOT intended to be used in place of a patient warming device. Follow patient warming device manufacturer's instructions for use to prevent patient hyperthermia or thermal insult. <p>Warnings:</p> <ul style="list-style-type: none"> The Patient Return Electrode utilizes current limiting technology that automatically limits the circuit current in the event the area of contact between the electrode and the patient drops below predetermined levels. Therefore, a CQM generator alarm will not sound if the patient is not in contact with the electrode. The current limiting technology built into the electrode is designed to prevent pad site burns by automatically limiting the circuit current in the event the area of contact between the patient and the electrode drops below predetermined levels, resulting in a noticeable reduction in the electrosurgical effect. Use the lowest possible power settings to achieve the desired effect. Use only with isolated generators. The Universal Plus or the Universal Dual Plus Patient Return Electrode may be used with patients weighing ≥ .35 kg (0.8 lb). There is no maximum weight restriction. Do not place the Patient Return Electrode directly onto a metal surface. Do not use the Patient Return Electrode with specialty modes on generators. For example, HIGH CUT or ENDO CUT modes when using the ERBE ICC 200, 300, 350, or VIO series of electrosurgical generator(s). Doing so may result in a different electrosurgical effect than intended. When the generator is operational, keep active accessories away from the patient and the Patient Return Electrode, when not in use, or store in an electrically isolated container. Do not use the Patient Return Electrode with any cuts or tears that expose the Polymer until properly repaired. The risk of igniting flammable gases or other materials is inherent in electrosurgery and cannot be eliminated by device design. Precautions must be taken to avoid contact of flammable materials and substances with electrosurgical electrodes, whether they are in the form of an anesthetic or skin preparation agent, are produced by natural processes within body cavities, or originate in surgical drapes, tracheal tubes, or other materials. To avoid the risk of tracheal fires, never use electrosurgery to enter the trachea during tracheostomy procedures. For surgical procedures where the high frequency current could flow through parts of the body having a relatively small cross sectional area (e.g. circumcisions), the use of bipolar techniques may be desirable in order to avoid unwanted tissue damage. Failure to cap the extra cord of the Universal Dual Plus Patient Return Electrode during use may result in an electrical hazard. <p>Cautions:</p> <ul style="list-style-type: none"> Refer to the generator manufacturer's operating manual for proper usage of electrosurgical equipment. As with other medical devices, avoid pooling of chemicals or fluids that may contribute to patient skin breakdown or pressure ulcer formation. Avoid placing the patient directly on the rigid corner or cable of the Patient Return Electrode to avoid possible pressure points. Follow proper nursing practices and procedures to protect the patient from pressure ulcers. Do not sterilize the Patient Return Electrode. Avoid sharp objects and edges at all times. Never use disinfectant containing more than 70% alcohol to clean the Patient Return Electrode. It may cause hardening of the outer skin. Only physicians familiar with electrosurgical effects on surgical patients should operate this device. 	<ul style="list-style-type: none"> The use and proper placement of a patient return electrode is a key element in the safe and effective use of monopolar electrosurgery, particularly in the prevention of burns. Follow directions and recommended practices for the preparation, placement, surveillance, and use of this patient return electrode. To reduce the chance of alternate site burns due to alternate current pathways, the patient, when practical should not be allowed to come into contact with metal parts that are grounded. It is recognized that this recommendation may not be practical during certain procedures; however, to maximize patient safety during the use of electrosurgical devices, such practices should be minimized. Electrosurgical cables should be positioned to minimize contact with the patient and avoid contact with other leads to avoid adversely influencing the operation of other electronic equipment. Apparent low power output or failure of the electrosurgical equipment to function correctly at normal settings may indicate faulty application of the patient return electrode or failure of an electrical lead. Do not increase power output before checking for obvious defects or misapplication. For monopolar surgery, effective contact between the patient and the patient return electrode must be verified whenever the patient is repositioned. Use of electrosurgical devices can cause severe electromagnetic interference in other devices, for patients with cardiac pacemakers, electrically conductive or other active implants, a possible hazard exists due to the concentration or re-direction of HF currents. The pacemaker or other active implant may be damaged due to the interference of HF currents. In case of doubt, approved qualified advice should be obtained from the device manufacturer. As with any patient return electrode, skin-to-skin contact (for example between the arms and body of the patient) should be avoided, for example by insertion of dry gauze. When high frequency surgical equipment and physiological monitoring equipment are used simultaneously on the same patient, any monitoring electrodes should be placed as far as possible from the surgical electrodes. Needle monitoring electrodes are not recommended. In all cases, monitoring systems incorporating high frequency current limiting devices are recommended. <p>Compatibility Statement:</p> <ul style="list-style-type: none"> Megadyne recommends the use of the Patient Return Electrode with the Megadyne electrosurgical generator and Megadyne accessory devices (i.e. EZ-CLEAN electrodes and pencils). <p>RF Voltage Rating For Neutral Electrodes Compliant with IEC 60601-2-2 requirements at frequency ratings from 300 to 1,000 kHz. RF Voltage Ratings ≤ 5.5 peak kV, consult electrosurgical generator specifications.</p> <p>U.S. Patents 6,053,910 6,214,000 6,454,764 6,544,258 6,582,424 6,666,859 D707,828, D709,205 Other U.S. and foreign patents issued and pending.</p> <p>The color green as used on Megadyne's E-Z Clean™ brand non-stick electrosurgical tips and brand patient return electrode is itself a registered trademark of Megadyne Medical Products Inc.</p>	
3						
4						

English

THIS DRAWING WAS RE-CREATED IN ADOBE INDESIGN

SIZE	PART NUMBER	REV.
B	3000314-01	003
DO NOT SCALE PRINT	SCALE N/A	SHEET 2 OF 27

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THIS DRAWING WAS RE-CREATED IN ADOBE INDESIGN

French EU

SIZE	PART NUMBER	REV.
B	3000314-01	003
DO NOT SCALE PRINT	SCALE N/A	SHEET 5 OF 27

1

2

3

4

**ETHICON
Megadyne™**

Megadyne Medical Products Inc.
11506 South State Street
Draper, Utah 84020 USA
Tel.: +1 (801) 576-9669

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