


AMT TLC[®] Self-Retaining Retractor System



KEY FEATURES:

- Low profile concave structure allows maximum access to operative site.
- Malleable hinges allow the retractor ring to adapt to the patient's anatomy quickly and easily.
- 6 different shapes and sizes to choose from.
- These retractors and stays are commonly used in colorectal, vaginal, penile, hand, and any other surface tissue surgeries.

Patents: 5,951,467, 6,090,043, 414,265

- Ribbed Elastic Stay Hooks fit snugly within the notched perimeter of the retractor frames. Stays range in size and bluntness.
- Finger grips on the distal portion of the stays make it easier for placement and repositioning.




Call one of our customer service representatives to place an order toll free at
1-800-869-7382
or email us at customerservice@appliedmedical.net

Applied Medical Technology, Inc., 8000 Katherine Blvd., Brecksville, OH 44141 USA

AMT TLC[®] Self-Retaining Retractor System

The TLC[®] Self-Retaining Retractor System has been designed so the surgeon can gently retract and secure the delicate soft tissue away from the operative surface. It allows the surgeon to perform the procedure without stay-sutures or assistance. This controlled versatility will give the surgical assistant the freedom to do other tasks, which will broaden the range of procedures the surgeon can perform single-handedly.



Simple and Versatile

- Affords the surgeon clear, unimpeded access to the surgical site.
- Offers increased control, support and exposure in multiple planes of retraction.
- Molded disposable plastic ring is cost-effective by reducing labor, cleaning and re-sterilization costs.




- Advanced one-piece construction eliminates the risk of detachment.
- Low-profile concave structure increases the accessibility of the stays for ease of placement.
- Ribbing on stays provides a secure fit into the ring, virtually eliminating slippage and risk of popping-out.
- Large finger grip on the stay enhances holding power, affording the surgeon increased control and safety.

Visit us on the web at www.appliedmedical.net


Applied Medical Technology, Inc., 8000 Katherine Blvd, Brecksville, Ohio 44141 USA

Specs - Dimensions in CM


Disposable retractor rings are pre-sterilized and packaged singly in Tyvek pouches for easy delivery to the sterile surgical field.




TLC200




TLC300




TLC400



TLC500



TLC600



TLC700

Call one of our customer service representatives to place an order toll free at
1-800-869-7382
or email us at customerservice@appliedmedical.net

AMT TLC® Self-Retaining Retractor System

Ribbed Elastic Stay Hooks

Sterile 4, 8 or 50 Packs

TLC150 Blunt Tip (12mm)

TLC110 Sharp Tip (5mm)

TLC100 Semi-Blunt Tip (5mm)

TLC090 Sharp Tip (3mm)



WILSON PENILE IMPLANTATION SYSTEM



Steven K. Wilson, M.D., F.A.C.S.,

Medical director of Southwest Impotency Center. Dr. Wilson is board certified in Urology and a clinical professor of Urology at the University of Arkansas. He has performed over 6,000 penile implant surgeries.

TLC5042

Wilson Penile Implantation System

Sterile Kit Contents:

- 6 - Blunt Tip (12mm) Stay Hooks
- 1 - Sharp Tip (5mm) Stay Hook
- 1 - Elevation Strap
- 1 - Retractor Frame (30.4cm x 19cm)

TLC5040

Wilson Scrotal Pak

Sterile Kit Contents:

- 6 - Blunt Tip (12mm) Stay Hooks
- 1 - Sharp Tip (5mm) Stay Hook
- 1 - Elevation Strap

Call one of our customer service representatives to place an order toll free at
1-800-869-7382

or email us at customerservice@appliedmedical.net

Visit us on the web at www.appliedmedical.net

Applied Medical Technology, 8000 Katherine Blvd, Brecksville, Ohio 44141 USA