

EASYFEEL

Indications:

The EASYFEEL kit is indicated for use in the fields of reconstructive or cosmetic surgery in order to increase the volume of intra and sub-dermal space deficiencies, in all parts of the body via the addition of adipose tissue.

The use of the kit is more particularly indicated in facial rejuvenation or reconstruction.

Materials:

The EASYFEEL kit comprises:

- 1 infiltration cannula, which can be used with 10 to 50cc luer-lock syringes
- 1 suction cannula **with specific small holes** (should be used with a 10cc luer-lock syringe)
- 4 caps for Becton Dickinson 10cc syringes, in order to stay **in closed circuit** during centrifugation
- 4 luer-lock stoppers to close the bottom of syringes during centrifugation
- 2 luer-lock adaptators for fat transfer to a smaller syringe
- 2 injection cannulas mainly dedicated for face lipofilling

CE Marking:

IDM 93/42/EEC "Council Directive 93/42/EEC of 14 June 1993 concerning medical devices »"
Classification : Medical Device (Annexe II) (BSI Cert. No. CE 564067)

Packaging:

The EASYFEEL kit is sterile and conditioned with double package and sold by units of 5.
Every kit contains all the elements for suction and injection **except the syringes**.

CHEAP, EFFICIENT and SAFE



- Infiltration cannula (x1)
- Multi holes sampling cannula
diam 1.7mm x 2.1mm 150mm (x1)
- Injection cannula
 - o 0.9mm x 1.1mm length 50mm
 - o 0.9mm x 1.1mm length 70mm
- Caps for BD 10cc syringes (x4)
- Luer Lock adaptators (x2)
- Luer Lock stoppers (x4)