



M.D. RESOURCE

medical device resource corporation



***Since 1985,
your source for
medical device
engineering,
manufacturing
and distribution.***

ISO 13485
CERTIFIED

This is just a sample of the many products we have available

Systems



LS2 and LS2DP (on cover)

Liposuction System, shown with MicroAire PAL Console. Has LED digital display, 2 rotary vane oilless motors, 29 in Hg, pneumatic footswitch in single (LS2) or dual port (LS2DP) systems. Both systems have IV pole, dual canister holder, dual footpedal capability, or optional wireless footpedal. Shown with TMD (tumescent measuring device). Dual port provides additional vacuum source if needed, as well as for quick canister changes.

System has quad filter system; IEC (power) outlets for your additional systems like the PAL and infiltration pump; as well as a wire bin for footpedals, cords, etc.



The old Flagship

LS1000, flagship of the rotary vane systems. Many still in the field, in service. You may recognize it as one of your early liposuction systems.

TMD

Tumescent Measuring Device

The convenient pole top mount delivers effortless measurement of total tumescent delivery. Weight lost from hanging tumescent fluid bag is digitally displayed with ± 1 ml accuracy. The very low current circuitry enables the TMD to be "piggy backed" to our K-Pump or BPI infiltrators or it can be used as a "standalone" device.



Infusion Pump

'K' Pump. This quiet and efficient pump can be mounted on the LS2, the PAL console, or on the NX3. Infiltrate from 0 to 1000cc/min; has dual footpedal capability; three roller peristaltic head; safety lever and additional power connector for TMD (Tumescent Measuring Device). Uses AFTD/AFTDY autofuse tubing.



— **contact us for a complete selection!**

Systems

TOUCH UP

Small, quiet, very portable piston pump aspirator for follow up procedures or fat transfer. 27 in Hg, pneumatic footpedal, vacuum control, valve canister ring, 20 lbs.



NX 3

Super quiet piston pump aspirator, 29 in Hg, dual pneumatic footswitch, single port, single control, 45 lbs., lightweight full powered aspirator.

NX3 shown here with mounted k-pump infiltration pump, dual canisters and TMD.



RM 2000

Quietest Aspirator on the market. 2 rotary vane oilless motors, 29 in Hg, Dual pneumatic footswitch, single port, single valve, dual canister and dual footpedal capability.



TITAN Rotary Vane

Whisper Silent. Low profile design of the RM 2000. Single port, single valve. Continuous switch and dual footswitch control.



Infiltrator BPI

Bi-phasic peristaltic infiltration pump for infusion, smaller and quieter than the k-pump. Infusion rates up to 900 ml/min. Comes with dual footpedal capability and works with the TMD as well. Uses SOFTUBE tubing equivalent to CG-Inf-T.

This is just a sample of the many products we have available

Accessories



AFTD and AFTDY Autofuse Tubing (Disposable)

Also available in
"Y" connection.



LS7010 and LS7012

Standard liposuction
tubing, very pliable not
stiff, 3/8" ID, 8' and 10'
lengths



SOFTUBE

Soft, pliable infiltration tubing for
use with the BPI Infiltration pump,
equivalent to CG Inf-T.



ITDY & ITD

Infiltration tubing
available in both
"Y" and normal,
equivalent to
Lamis tubing

LS7020 and LS7030

Implosion proof canisters, 2000cc
and 3000cc. No inserts, direct
and easy view and measurement.
Disposable.



FL3000

Absolute Biofilter, industry
standard 3/8" connector, fits all
liposuction aspirators.

