


Height (mm)	Diameter (mm)	Volume (ml)	Reference
25	21	6	LS 01 025
34	27	12	LS 01 034
40	30	17	LS 01 040
50	33	27	LS 01 050
53	38	38	LS 01 053

Not every device will conform to the dimensions given. Slight variations may occur. On demand, production of a textured surface (reference LST).

SHAPE AND CONSISTENCY OF A TESTICLE

Various sizes corresponding to patient morphology



Our sterile and single use testicular implants are sold by unit.

CONTINUITY Manufacturing process developed throughout the last 26 years.

RELIABILITY Chemical and biological evaluation of raw materials then microbiological evaluation of each batch. Visual and technical controls on 100% of implants at each step of the manufacturing process.

TRACEABILITY Each implant identification number inseparable from its reference allows traceability from the raw materials up to the implantation. We introduce our new bar code system GSI-128 to improve everyday life.

SERENITY Highly qualified staff working in clean room.



CE 0483 in accordance with
Medical Device Directive 93/42 EEC



Rubanksgatan 8
SE-741 71 Knivsta
SWEDEN

Tel: +46(0)18-54 54 00
Fax: +46(0)18-51 50 50
E-mail: promeduc@promeduc.se
Web: www.promeduc.se