


HUNTER™ ACTIVE AND PASSIVE TENDON IMPLANTS


*a consistently
reliable method
of salvaging
scarred tendon
systems*


- Designed for temporary use in reconstruction of extension and flexor tendons based on a two stage procedure.
- Rods create a pseudosynovial sheath which later lubricates and nourishes an autogenous graft.
- Manufactured with polyester reinforcement providing added strength and reducing risk of detachment and migration.
- Bone screw plate available providing firm fixation of distal attachment further reducing the risk of detachment and migration
- Active tendons allow use of hand during rehabilitation thereby:
 - increasing patient motivation to comply with therapy
 - reducing loss of range of motion before autogenous graft implantation
 - reducing proximal juncture adhesions at the site of motor tendon sheath
 - causing softening of the connective tissue and joints
 - allowing active evaluation of the pulley system before autogeneous graft placement





PROMEDUC MEDICAL AB
Börjég. 50, 752 29 UPPSALA
☎ 018 - 54 54 00
e-mail: promeduc@promeduc.se
web: www.promeduc.se


HUNTER™ ACTIVE AND PASSIVE TENDON IMPLANTS


<i>HUNTER™ Active Tendon Implant (4 mm Width)</i>		<i>Order</i>
16 cm Long		AT164000
18 cm Long		AT184000
20 cm Long		AT204000
22 cm Long		AT224000

<i>HUNTER™ Active Tendon Implant PC (4mm Width)</i>		
16 cm Long		ATPC1640
18 cm Long		ATPC1840
20 cm Long		ATPC2040
22 cm Long		ATPC2240
26 cm Long		ATPC2640

<i>HUNTER™ Active Tendon Implant DC (4 mm Width)</i>		
16 cm Long		ATDC1640
18 cm Long		ATDC1840
20 cm Long		ATDC2040
22 cm Long		ATDC2240
26 cm Long		ATDC2640

<i>HUNTER™ Active Tendon Implant BC (4 mm Width)</i>		
16 cm Long		ATBC1640
18 cm Long		ATBC1840
20 cm Long		ATBC2040
22 cm Long		ATBC2240
26 cm Long		ATBC2640

<i>HUNTER™ Active Tendon Implant (25 cm Length)</i>		
3 mm Wide		PT300000
4 mm Wide		PT400000
5 mm Wide		PT500000
6 mm Wide		PT600000

<i>HUNTER™ Tendon Rod (24.5 cm Length)</i>		
2 mm Wide		TR200000
3 mm Wide		TR300000
4 mm Wide		TR400000
5 mm Wide		TR500000
6 mm Wide		TR600000



A Wright Medical Group Company

5677 AIRLINE ROAD • ARLINGTON, TN, USA 38002
800-238-7117 • 901-867-9971 • WWW.WMT.COM

©2001 Wright Medical Technology, Inc.
All rights reserved.

HUNTER™ is a trademark of
Wright Medical Technology, Inc.

S0773-798
(Rev. 3-01)