

Restoring a body.  Restoring a life.



Shape an
ideal path to
reconstruction.

MENTOR® CPX™4
Breast Tissue
Expander

NOT FOR DISTRIBUTION IN USA.
APPROVED FOR EMEA REGION.

 **MENTOR®**
Make life more beautiful™



CHOOSE
THE TISSUE
EXPANDER
THAT
CONSIDERS
YOUR NEEDS.

Unique design features
and enhancements
optimize performance
for expansion results
you and your patients
will welcome.

**MENTOR® CPX™4 BREAST
TISSUE EXPANDERS**

EXPERIENCE THE DIFFERENCE THROUGHOUT THE EXPANSION CONTINUUM.

INTRAOPERATIVE

*Smooth, flat injection dome and enhanced shell pliability**

- Helps ease insertion and removal through incision

Optional suture tabs at 3, 6 and 9 o'clock positions

- Stabilize placement to help predict and control the location of the reconstructed inframammary fold
- Allow left to right alignment in bilateral procedures

DURING THE EXPANSION

BufferZone™ Self-Sealing Patch and posterior Dacron® Patch

- Focus expansion towards the lower pole of the breast

More pliant proprietary BufferZone™ Patch and smooth injection dome surface

- Designed for enhanced patient comfort
- Prevent leakage in the event of an accidental needle puncture

*48% stronger magnet**

- Designed to indicate the optimal position of the fill port through the CENTERSCOPE™ Magnetic Detection Device

RELIABLE 2-STAGE RECONSTRUCTION

BufferZone™ Self-Sealing Patch and Dacron® Patch

- Drive lower pole expansion to mimic the shape of the natural breast

SILTEX® Texture

- Prevents pooling or sequestration of body fluids
- Provides disruptive surface for collagen interface

Implant Matching

- Matching range of definitive implants to be used after tissue expansion


*Compared with the MENTOR® CPX2™ and CPX3™ Breast Tissue Expanders.

MENTOR® CPX™4 Breast Tissue Expanders


Optionally available with suture tabs at 3, 6 and 9 o'clock positions

FOR MORE INFORMATION OR
TO PLACE AN ORDER,
visit www.mentorwwllc.eu


THREE CONTOURED SHAPES ACCOMMODATE ALL YOUR NEEDS.




MENTOR® CPX™4 Breast Tissue Expander Style 8100 - Low Height					
CATALOG	VOLUME	WIDTH	HEIGHT	PROJECTION	
3548121	250 cc	11.4 cm	8.1 cm	6.1 cm	
3548122	350 cc	12.7 cm	9.4 cm	6.5 cm	
3548123	450 cc	14.0 cm	10.2 cm	7.1 cm	
3548124	550 cc	15.0 cm	10.9 cm	7.4 cm	
3548125	650 cc	15.7 cm	11.2 cm	7.9 cm	
3548126	750 cc	16.5 cm	11.9 cm	8.1 cm	




MENTOR® CPX™4 with Suture Tabs Breast Tissue Expander Style 9100 - Low Height					
CATALOG	VOLUME	WIDTH	HEIGHT	PROJECTION	
3549121	250 cc	11.4 cm	8.1 cm	6.1 cm	
3549122	350 cc	12.7 cm	9.4 cm	6.5 cm	
3549123	450 cc	14.0 cm	10.2 cm	7.1 cm	
3549124	550 cc	15.0 cm	10.9 cm	7.4 cm	
3549125	650 cc	15.7 cm	11.2 cm	7.9 cm	
3549126	750 cc	16.5 cm	11.9 cm	8.1 cm	




MENTOR® CPX™4 Breast Tissue Expander Style 8200 - Medium Height					
CATALOG	VOLUME	WIDTH	HEIGHT	PROJECTION	
3548221	275 cc	10.7 cm	9.3 cm	6.2 cm	
3548222	350 cc	11.7 cm	10.0 cm	6.6 cm	
3548223	450 cc	12.7 cm	10.8 cm	7.0 cm	
3548224	550 cc	13.5 cm	11.7 cm	7.4 cm	
3548225	650 cc	14.6 cm	12.6 cm	7.6 cm	
3548226	800 cc	15.6 cm	13.3 cm	8.0 cm	



MENTOR® CPX™4 with Suture Tabs Breast Tissue Expander Style 9200 - Medium Height					
CATALOG	VOLUME	WIDTH	HEIGHT	PROJECTION	
3549221	275 cc	10.7 cm	9.3 cm	6.2 cm	
3549222	350 cc	11.7 cm	10.0 cm	6.6 cm	
3549223	450 cc	12.7 cm	10.8 cm	7.0 cm	
3549224	550 cc	13.5 cm	11.7 cm	7.4 cm	
3549225	650 cc	14.6 cm	12.6 cm	7.6 cm	
3549226	800 cc	15.6 cm	13.3 cm	8.0 cm	



MENTOR® CPX™4 Breast Tissue Expander Style 8300 - Tall Height					
CATALOG	VOLUME	WIDTH	HEIGHT	PROJECTION	
3548321	250 cc	10.1 cm	10.7 cm	5.6 cm	
3548322	350 cc	11.3 cm	11.8 cm	6.0 cm	
3548323	450 cc	12.3 cm	12.9 cm	6.5 cm	
3548324	550 cc	13.2 cm	13.8 cm	6.9 cm	
3548325	650 cc	14.0 cm	14.6 cm	7.3 cm	
3548326	750 cc	14.6 cm	15.3 cm	7.6 cm	
3548327	850 cc	15.4 cm	15.9 cm	7.9 cm	



MENTOR® CPX™4 with Suture Tabs Breast Tissue Expander Style 9300 - Tall Height					
CATALOG	VOLUME	WIDTH	HEIGHT	PROJECTION	
3549321	250 cc	10.1 cm	10.7 cm	5.6 cm	
3549322	350 cc	11.3 cm	11.8 cm	6.0 cm	
3549323	450 cc	12.3 cm	12.9 cm	6.5 cm	
3549324	550 cc	13.2 cm	13.8 cm	6.9 cm	
3549325	650 cc	14.0 cm	14.6 cm	7.3 cm	
3549326	750 cc	14.6 cm	15.3 cm	7.6 cm	
3549327	850 cc	15.4 cm	15.9 cm	7.9 cm	

Individual implant dimensions may vary slightly in products of this type. Not all units conform exactly to the dimensions noted above.

Important Safety Information:

MENTOR® CONTOUR PROFILE™ Tissue Expanders and MENTOR® Smooth Round Tissue Expanders with Remote Dome are used for breast reconstruction following mastectomy, correction of an underdeveloped breast, scar revision and tissue defect procedures. These expanders are intended for temporary subcutaneous or submuscular implantation; they should be used within a time frame determined by the physician to achieve the clinically desired degree of tissue expansion. CONTOUR PROFILE™ Tissue Expanders are devices that contain magnetic injection domes and are NOT MRI compatible. Do not use the CONTOUR PROFILE™ Tissue Expander in patients where an MRI may be needed. DO NOT use the CONTOUR PROFILE™ Tissue Expander in patients that have a previously implanted device that could be affected by a magnetic field. The device could be moved by the MRI causing pain or displacement, potentially resulting in a revision surgery. The incidence of extrusion of the expander has been shown to increase when the expander has been placed in injured areas: scarred, heavily irradiated or burned tissue, crushed bone areas or where severe surgical reduction of the area has previously been performed. Your patient needs to be informed and understand the risks and benefits of MENTOR® Tissue Expanders, and provided with an opportunity to consult with you prior to deciding on surgery. For detailed indications, contraindications, warning and precautions associated with the use of all MENTOR® Implantable Devices, please refer to the Product Insert Data Sheet provided with each product, or review the Important Safety Information provided www.mentorwwllc.eu.

Operations Europe, Middle East, Africa

Germany, Austria & Switzerland
Mentor Deutschland GmbH
Tel: +49 811 600 50 0
Fax: +49 811 600 50 11

South Africa
Johnson & Johnson Medical (Pty) Ltd
Tel: +27 11 265 1135
Fax: +27 11 265 1114

**Croatia, Macedonia, Serbia,
Slovenia & Bulgaria**
Johnson & Johnson S.E. d.o.o.
Tel: +386 1 40 11 800
Fax: +386 1 40 11 801

Spain & Portugal
Mentor Medical Systems Iberica, S.L.
Tel: +34 91 722 8000
Fax: +34 91 728 3367

Turkey
Johnson & Johnson Medical Turkey
Tel: +90 216 538 21 21
Fax: +90 216 538 21 30

**Czech Republic, Slovak Republic,
Hungary & Romania**
Johnson & Johnson s.r.o.
Tel: +420 227 012 371
Fax: +420 227 012 300

Italy, Greece & Israel
Johnson & Johnson Medical S.p.a.
Tel: +39 06.91194500
Fax: +39 06.91194505

Middle East, Egypt & Pakistan
Johnson & Johnson Middle East Inc
Tel: +971 4 429 72 00
Fax: +971 4 299 43 30

**Poland, Lithuania,
Latvia & Estonia**
Johnson & Johnson Poland sp. z o.o.
Tel: +48 22 237 83 00
Fax: +48 22 237 83 53

France & Benelux
Perouse Plastique, SAS
Tel: +33 1 55 00 22 45
Fax: +33 1 55 00 28 11

Maghreb & Sub Saharan Africa
Johnson & Johnson Professional Export
Tel: +44 (0)1506 594 500
Fax: +44 (0)1506 594 731

**Sweden, Denmark, Norway,
Finland & Iceland**
Johnson & Johnson AB
Tel: +46 8 626 2200
Fax: +46 8 754 5850

United Kingdom & Ireland
Johnson & Johnson Medical Ltd
Tel: +44 (0)800 973 082
Fax: +44 (0)134 486 4122

Russia & CIS
Johnson & Johnson LLC
Tel: +7 495 580 7777
Fax: +7 495 580 7878

Corporate Headquarters
Santa Barbara, CA 93111 USA
Tel: +1 805 879 6000
www.mentorwwllc.com