


BIOVER

DISPOSABLE VESSEL CLAMP



BIOVER® AG
www.biover.ch

Distributed by:

PROMEDUC
surgical

Promeduc Surgical AB
Rubanksgatan 8, 741 71 Knivsta, SWEDEN
tel +46-18-54 54 00
promeduc@promeduc.se www.promeduc.se

DISPOSABLE MICROVASCULAR CLAMPS

The BIOVER[®] microvascular clamp is a high precision instrument and has been developed under the proper guidance of many international well-known surgeons.



- Highest quality (Polycarbonate and stainless steel)
- Single use only (to be discarded after use)
- Guaranteed sterile
- Cheap in cost
- Three different sizes are available for arteries and four sizes for veins:

- for arteries:


- Ø less than 1 mm TKS
- Ø 1 to 2 mm TKM
- Ø 2 to 4 mm TKL

- for veins:

- Ø 0.2 to 1 mm TKFV
- Ø less than 1 mm TKS
- Ø 1 to 2 mm TKMV
- Ø 2 mm TKLV



Actual Size



Ordering Information

Size/Cat.-No.	A	B	C	D	Weight	Clipping power
Arteries (yellow clamp)						
TKS-1	10.0 mm	4.5 mm	2.5 mm		0.06 g	40 g
TKS-2				~ 7.5 mm	0.13 g	
TKM-1	11.5 mm	5.5 mm	3.0 mm		0.12 g	60 g
TKM-2				~12.5 mm	0.25 g	
TKL-1	15.0 mm	6.5 mm	3.5 mm		0.25 g	120 g
TKL-2				~15.0 mm	0.50 g	
Veins (green clamps)						
TKFV-1	7.0 mm	4.0 mm	2.3 mm		0.03 g	15 g
TKFV-2				~8.0 mm	0.06 g	
TKSV-1	10.0 mm	4.5 mm	2.5 mm		0.06 g	20 g
TKSV-2				~7.5 mm	0.13 g	
TKMV-1	11.5 mm	5.5 mm	3.0 mm		0.12 g	30 g
TKMV-2				~12.5 mm	0.25 g	
TKLV-1	15.0 mm	6.5 mm	3.5 mm		0.25 g	40 g
TKLV-2				~15.0 mm	0.50 g	

Conditions:

- Single or double clamps, sterile, ready for use.
- Delivery: 10 single or 10 double clamps per box.
- Each box easily identifiable by colour (each type a different colour).

Warning:

Do not use the BIOVER vascular clamps for the central blood circulation and not for the central nervous system.