

# plasmasure

OPTIMIZED SOLUTION FOR PLATELET-RICH PLASMA PREPARATION

plasmasure kit is a sterile single use kit for the preparation of autologous Platelet-Rich Plasma solution.

## ADVANTAGES

A complete blood treatment kit including **all the elements** required to perform a PRP preparation and injection

Blood remains in a **closed circuit** device, reducing the risk of infection for patients

**User friendly** device

A pure PRP containing more than **80%** of platelets

Highest Platelet Capture Efficiency within a **single centrifugation**

**20ml** of sampled blood = **5 to 7ml** of injectable PRP

## STEMCIS

A BRAND OF DMS BIOTECH

Distributed by:

**PROMEDUC**  
*surgical*

Promeduc Surgical AB | Rubanksgatan 8, 741 71 Knivsta, SWEDEN  
tel +46-18-54 54 00 | promeduc@promeduc.se | www.promeduc.se

# A SINGLE USE STERILE KIT

**plasmasure** kit contains all the necessary items from blood collection to PRP injection.



- 1 syringe (20ml)
- 1 syringe (10ml)
- 1 syringe (5ml)
- 1 syringe (1ml)
- 1 mesotherapy needle 30G L4mm
- 1 mesotherapy needle 27G L13mm
- 1 mesotherapy needle 30G L13mm
- 1 mesotherapy needle 27G L6mm
- 1 needle 21G L50mm
- 1 needle 22G L50mm

- 2 transfer cannulas (C)
- 1 butterfly needle 21G L19mm
- 1 centrifuge tube (T)
- 1 balance tube
- 2 connectors
- 1 nut for centrifuge rotor
- 2 Luer Lock stoppers
- 1 Sodium Citrate 3.8% anticoagulant solution
- 1 sterile surgical drape

