



Smooth button  
Minimizes risk of  
damaging gloves

Smoke evacuation

RetraLight® is a single  
use retractor with  
integrated LED light  
and smoke evacuation.


Antislip design  
Reduces fatigue during lift

RetraLight® is  
designed specifically  
for breast surgery  
together with leading  
plastic surgeons.

Easy accessibility  
Quick recycling of the battery

Ergonomic handgrip  
Fits all sizes of hands comfortably

With or without teeth  
Available in two versions



Smooth bottom section  
For implant and tissue safety

Cordless illumination  
Using a powerful LED

RL4000 RetraLight® LED Retractor  
(with teeth)  
RL4001 RetraLight® LED Retractor  
(without teeth)

**PROMEDUC**  
*surgical*

Promeduc Surgical AB, Rubanksgatan 8 SE-741 71 Knivsta | [www.promeduc.se](http://www.promeduc.se)