



AESTHET/C
GROUP

Fat transfer packs

Sterile components
&
single use

order@promeduc.se

018-54 54 00

AESTHET/C

NANOGRAFT - Volume 1-10mL

Conçu pour une procédure qui requiert un **très petit volume** de tissu adipeux

*Designed for a procedure that requires a **very small volume** of adipose tissue*



REF. PACK-NALT

Description	REF.	Pc(s)
Canule d’infiltration Ø2.10mm x 150mm – 10 orifices <i>Infiltration cannula</i>	CI21150	1
Canule de prélèvement Ø2.10mm x 150mm – 8 orifices <i>Harvesting cannula</i>	CP141508H	1
Micro-canule d’injection 22G Ø0.70mm x 50mm – 1 orifice <i>Injection cannula</i>	MI2250KR	1
Micro-canule d’injection 22G Ø0.70mm x 70mm – 1 orifice <i>Injection micro cannula</i>	MI2270KR	1
Raccord de transfert FLL/FLL Øintern. 2.4mm <i>Connector</i>	FD0000000-LLF24	1
Raccord de transfert FLL/FLL Øintern. 1.4mm <i>Connector</i>	FD0000000-LLF14	1
Raccord de transfert FLL/FLL Øintern. 1.2mm <i>Connector</i>	FD0000000-LLF12	1
Seringues 10ml LL BD® <i>Syringes</i>	300912	2
Seringues 1ml LL BD® <i>Syringes</i>	309628	2
Bouchons LL <i>Caps</i>	FD0001711-MC	2
Aiguille 14G x40mm <i>Needle</i>	8300013707 (AF1440)	1
Aiguille 21g x40mm <i>Needle</i>	4710008040	1

MEDIUMGRAFT - Volume 10-30mL

Conçu pour une procédure qui requiert un **volume moyen** de tissu adipeux

*Designed for a procedure that requires a **medium volume** of adipose tissue*



REF. PACK-MELT

Description	REF.	Pc(s)
Canule d’infiltration Ø2.10mm x 150mm – 10 orifices <i>Infiltration cannula</i>	CI21150	1
Canule de prélèvement Ø2.10mm x 150mm – 6 orifices spiral <i>Harvesting cannula</i>	CP141506HS	1
Canule d’injection 14G Ø2.10mm x 150mm – 1 orifice 90° <i>Injection cannula</i>	CF201590	1
Canule d’injection 17G Ø1.45mm x 90mm – 1 orifice 90° <i>Injection cannula</i>	CF14790	1
Raccord de transfert FLL/FLL Øintern. 1.2mm <i>Connector</i>	FD0000000-LLF12	2
Raccord de transfert FLL/FLL Øintern. 1.4mm <i>Connector</i>	FD0000000-LLF14	2
Raccord de transfert FLL/FLL Øintern. 2.4mm <i>Connector</i>	FD0000000-LLF24	2
Seringues d’aspiration de 30ml avec piston autobloquant <i>Vacuum syringes with self-locking plunger</i>	ZZSQ30	4
Seringue 10ml LL BD® <i>Syringes</i>	305959	6
Seringues 1ml LL BD® <i>Syringes</i>	309628	6
Bouchons LL <i>Caps</i>	FD0001711-MC	4
Aiguille 14G x40mm <i>Needle</i>	8300013707 (AF1440)	2

Les tout-en-un pour les différentes étapes d'infiltration, prélèvement et de réinjection !
The all-in-one for the different step of infiltration, harvesting and reinjection !

MICROGRAFT - Volume 1-10mL

Conçu pour une procédure qui requiert un **petit volume** de tissu adipeux

*Designed for a procedure that requires a **small volume** of adipose tissue*



REF. PACK-MILT

Description	REF.	Pc(s)
Canule d’infiltration Ø2.10mm x 150mm – 10 orifices <i>Infiltration cannula</i>	CI21150	1
Canule de prélèvement Ø2.10mm x 150mm – 8 orifices <i>Harvesting cannula</i>	CP141508H	1
Micro-canule d’injection 18G Ø1.25mm x 90mm – 1 orifice <i>Injection micro cannula</i>	MI1890KR	1
Micro-canule d’injection 21G Ø0.81mm x 40mm – 1 orifice <i>Injection micro cannula</i>	MI2140KR	1
Micro-canule d’injection 21G Ø0.81mm x 70mm – 1 orifice <i>Injection micro cannula</i>	MI2170KR	1
Raccord de transfert FLL/FLL Øintern. 1.2mm <i>Connector</i>	FD0000000-LLF12	1
Raccord de transfert FLL/FLL Øintern. 1.4mm <i>Connector</i>	FD0000000-LLF14	1
Raccord de transfert FLL/FLL Øintern. 2.4mm <i>Connector</i>	FD0000000-LLF24	1
Seringue 10ml LL BD® <i>Syringes</i>	300912	4
Seringues 1ml LL BD® <i>Syringes</i>	309628	4
Bouchons LL <i>Caps</i>	FD0001711-MC	4
Aiguille 14G x40mm <i>Needle</i>	8300013707 (AF1440)	1
Aiguille 21g x40mm <i>Needle</i>	4710008040	1

MACROGRAFT - Volume 10-60mL

Conçu pour une procédure qui requiert un **volume important** de tissu adipeux

*Designed for a procedure that requires a **large volume** of adipose tissue*



REF. PACK-MALT

Description	REF.	Pc(s)
Canule d’infiltration Ø2.10mm x 150mm – 10 orifices <i>Infiltration cannula</i>	CI21150	1
Canule de prélèvement Ø3.0mm x 230mm – 24 orifices <i>Harvesting cannula</i>	CP302324H	1
Canule d’injection 14G Ø2.10mm x 150mm – 1 orifice 90° <i>Injection cannula</i>	CF201590	1
Canule d’injection 14G Ø2.10mm x 220mm – 1 orifice 90° <i>Injection cannula</i>	CF202290	1
Raccord de transfert FLL/FLL Øintern. 1.2mm <i>Connector</i>	FD0000000-LLF12	2
Raccord de transfert FLL/FLL Øintern. 1.4mm <i>Connector</i>	FD0000000-LLF14	2
Raccord de transfert FLL/FLL Øintern. 2.4mm <i>Connector</i>	FD0000000-LLF24	2
Seringues d’aspiration de 30ml avec piston autobloquant <i>Vacuum syringes with self-locking plunger</i>	ZZSQ30	6
Seringue 10ml LL BD® <i>Syringes</i>	305959	6
Bouchons LL <i>Caps</i>	FD0001711-MC	4
Aiguille 14G x40mm <i>Needle</i>	8300013707 (AF1440)	1